



# Ausrüstungen für medizinische Gase



## SOCKET LEAF - CZ (CSN 85 2762)



Organismo Notificato n. 0476  
Notified Body nr. 0476

### Produktbeschreibung

The socket of the new LEAF terminal unit line meets the needs of functionality, safety, innovation and design, in compliance with the CSN 85 2762 standard; The socket has a brass body and a male coupling locking system with steel balls; the coupling flange with the base is rotatable so you can position the terminal unit in the desired orientation; fixing to the base block is simple and maintenance is quick thanks to the replacement of a single pre-assembled cartridge; thanks to the built-in NFC system, traceability and maintenance management are facilitated via the MyBehringer App, which can be installed on

your smartphone.

### Normen

UNI EN ISO 9170-1 | UNI CEI EN ISO 13485 | UNI EN ISO 15001 | CSN 85 2762

### Komponenten

- Chemically nickel-plated brass body.
- Release ring.
- Stainless steel springs, balls, and filter.
- EPDM perox sealing O-ring.

### Wartungskit

CODE	GAS	DESCRIPTION
HK210CZ00	ALL	CARTUCCIA MANUT. PRESA LEAF CSN 85 2762 - N°11

### Technische Daten

CODE	GAS	FLOW RATE	MAXIMUM PRESSURE DROP	WEIGHT
HP301CZA4BEH	Air	40 l/min	0,15 bar	0,12 kg
HP301CZCDBEH	CO <sub>2</sub>	40 l/min	0,15 bar	0,12 kg
HP301CZNOBEH	N <sub>2</sub> O	40 l/min	0,15 bar	0,12 kg
HP301CZOXBEB	O <sub>2</sub>	40 l/min	0,15 bar	0,12 kg

## Technische Daten

CODE	GAS	FLOW RATE	MAXIMUM PRESSURE DROP	WEIGHT
HP301CZVABEH	Vac	25 l/min	0,15 bar	0,12 kg
HP301CZA8BEH	Air-800	350 l/min	0,70 bar	0,12 kg
HP301CZO9BEH	O2 93	40 l/min	0,15 bar	0,12 kg

## Technische bilder

