



Medical gas devices



Central ramp with non-return valve


 Organismo Notificato n. 0476
 Notified Body nr. 0476

Product description

Cylinder connection ramp made with a single brass manifold, with two or three places; specific gas connection; non-return valve built-in in the inlet connection. It is made of materials compatible with the gases and applicable standards in the medicinal gas sector. The valve bodies are made of brass. The identification

takes place through the indication of the name and/or symbol of the supplied gas and the production batch stamped on the manifold. The inlet and outlet connectors are made of brass with gas specific threads as per table. The wall fixing bracket is in galvanized steel.

Normatives

UNI EN ISO 7396-1 | UNI EN ISO 21969 | UNI EN ISO 15001

Components

- Body in brass.
- Inlet connection with non return valve in brass.
- Inox plate.


Maintenance kit

CODE	GAS	DESCRIPTION
HV039N02	O ₂	EQ. NON-RETURN VALVE W21,7x1/14" 'A'
HV039N09	N ₂ O	EQ. NON-RETURN VALVE G3/8" 'A'
HV039N05	N ₂	EQ. NON-RETURN VALVE W21,7x1/14" 'B'
HV039N10	CO ₂	EQ. NON-RETURN VALVE W27x2 'A'
HV039N06	Air	EQ. NON-RETURN VALVE W30x1/14" 'A'

Technical data

CODE	GAS	N. FOLDS	CONNECTION	WEIGHT	P1.max
HS0900X	O ₂	2	UNI-2	2,70 kg	200 bar
HS090NO	N ₂ O	2	UNI-9-DIN-11	2,70 kg	200 bar
HS090NI	N ₂	2	UNI-5	2,70 kg	200 bar
HS090CD	CO ₂	2	UNI-10	2,70 kg	200 bar
HS090AI	Air	2	UNI-6	2,70 kg	200 bar
HS0910X	O ₂	3	UNI-2	4,20 kg	200 bar
HS091NO	N ₂ O	3	UNI-9-DIN-11	4,20 kg	200 bar
HS091NI	N ₂	3	UNI-5	4,20 kg	200 bar
HS091CD	CO ₂	3	UNI-10	4,20 kg	200 bar
HS091AI	Air	3	UNI-6	4,20 kg	200 bar

Technical Images



Behringer S.r.l. via Gualco 4, 16165 Genova ITALY

T +39 010 8309103

F +39 010 8309183

e-mail info@behringer.it

C.F. e P.I. 01167840998

